



OCEAN King India

Home of Quality Products

Ocean FDM-M1 Biological 3D Printer



M1

Biological 3D printer

Forming technique	FDM
Print dimension:X	250mm
Print dimension:Y	150mm
Print dimension:Z	140mm
layer thickness(Z)	0.5mm
Precision	± 0.1mm/100mm
Speed	25~250mm/s
Power supply	VAC 50/60Hz 110V~220V
Output	120W
Weight	18kg
Overall dimension	(L)486mmX(W)377mmx(H)425mm
Material	Biological liquid materials

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471

Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co