



DJ-A Series Precise Electronic Balance

Balances with 0.1g resolution

Model	Capacity Division	Pan Size (mm)	Display Method	Power
DJ301A	300g/0.1g	120	LCD with backlight	AC/DC 6V/2.8Ah DC9V 500 mA
DJ501A	500g/0.1g			
DJ601A	600g/0.1g			
DJ1001A	1000g/0.1g			
DJ1201A	1200g/0.1g			

Balances with 0.01g resolution

Model	Capacity Division	Pan Size (mm)	Display Method	Power
DJ102A	100g/0.01g	120 weighing wind shield include	LCD with backlight	AC/DC 6V/2.8Ah DC9V 500 mA
DJ202A	200g/0.01g			
DJ302A	300g/0.01g			
DJ502A	500g/0.01g			
DJ602A	600g/0.01g			
DJ1002A	1000g/0.01g			
DJ2002A	2000g/0.01g			

Balances with 2000g/0.1g resolution

Model	Capacity Division	Pan Size (mm)	Display Method	Power
DJ2001A	2000g/0.1g	160*170	LCD with backlight	AC/DC 6V/2.8Ah DC9V 500 mA
DJ3001A	3000g/0.1g			
DJ5001A	5000g/0.1g			
DJ6001A	6000g/0.1g			

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471
Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co



Balances with 0.001g resolution

Model	Capacity Division	Pan Size (mm)	Display Method	Power
DJ103A	100g/0.001g	120 weighing wind shield include	LCD with backlight	AC/DC 6V/2.8Ah DC9V 500 mA
DJ103A-2	100g/0.002g			
DJ203A-2	200g/0.002g			
DJ203A-5	100g/0.005g			
DJ303A-5	300g/0.005g			

Common Features:

- * Large LCD display with auto - backlight
- * Selectable measure units: g.oz.ct.lb.tl.dwt
- * Auto-calibration through key pad
- * Auto-tare Auto-counting
- * Overload protection
- * Low voltage protection
- * AC/DC dual powered, Rechargeable battery powered for back up
- * RS-232 output
- * Extra display available
- * Scale dimension: 270*170*80 (mm)
- * Carton dimension: 600*380*510 (mm)/(6PCS/CTN)
- * Carton (N.W/G.W): 18kg/20kg

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471
Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co