



High quality VFD Display DJV

Balance with 0.1g resolution

Model	Capacity Division	Pan Size (mm)	Display Method	Power
DJ101V	100g/0.1g	120	VFD Display	AC220v/50Hz
DJ201V	200g/0.1g			
DJ301V	300g/0.1g			
DJ501V	500g/0.1g			
DJ601V	600g/0.1g			
DJ1001V	1000g/0.1g			
DJ1201V	1200g/0.1g			

Balance with 0.01g resolution

Model	Capacity Division	Pan Size (mm)	Display Method	Power
DJ102V	100g/0.01g	120weighing wind shield include	VFD Display	AC220v/50Hz
DJ202V	200g/0.01g			
DJ302V	300g/0.01g			
DJ502V	500g/0.01g			
DJ602V	600g/0.01g			
DJ1002V	1000g/0.01g			

Common Features:

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471

Email: info@oceankingindia.com, amitsood24@gmail.com

website: www.oceankingindia.com, www.goldtester.co



OCEAN King India

Home of Quality Products

- * VFD Display
- * Auto-calibration through key pad
- * Auto-tare Auto-counting
- * Overload protection
- * Selectable measure units:g.oz.kg.lb.tl.dwt
- * low voltage protection
- * RS-232 output
- * Scale dimension: 250*172*80 (mm)
- * Carton dimension: 600*360*300 (mm)/4(PCS/CTN)
- * Carton (N.W/G.W): 10.5g/11.5kg

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471

Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co