



OCEAN King India

Home of Quality Products



Silver Series DJKD

DJ-KD Series Precise Electronic Balance

Model	Capacity Division	Pan Size (mm)	Display Method	Power
DJ3KD	3kg/0.2g	260*210	LCD Display	6V/4AhDC9V 500mA
DJ6KD	6kg/0.5g			
DJ15KD	15kg/1g			
DJ30KD	30g/2g			

Common Features:

- * Fast and stable Display
- * Auto-calibration
- * Auto-tare
- * Counting Function
- * Selectable measure units:g.oz.kg.lb.tl.dwt
- * low voltage protection
- * AC/DC dual powered Rechargeable battery powered for back up
- * Daul span higher precision
- * Overload protection
- * Scale dimension: 310*265*88 (mm)
- * Carton dimension: 720*360*440 (mm)/4(PCS/CTN)
- * Carton (N.W/G.W): 24kg/26kg

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471
Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co