



OCEAN King India

Home of Quality Products



Single Display DJB

Balances with 0.1g resolution

Model	Capacity Division	Pan Size (mm)	Display Method	Power
DJ102B	100g/0.01g	115 weighing wind shield	LCD with backlight	AC/DC 6V/2.8Ah DC9V 300 mA
DJ202B	200g/0.01g			
DJ301B	300g/0.01g			
DJ501B	500g/0.01g			
DJ601B	600g/0.01g			
DJ1002B	1000g/0.01g			

Common Features:

- * Large LCD display with auto - backlight
- * Auto-calibration
- * Auto-counting
- * Auto-tare
- * Selectable measure units: g.oz.ct.lb.tl.dwt
- * Overload protection
- * Low voltage protection
- * RS-232 output
- * Extra display available
- * Scale dimension: 230*175*75 (mm)
- * Carton dimension: 600*380*510 (mm)/(6PCS/CTN)
- * Carton (N.W/G.W): 13.5kg/15kg

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471

Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co