



# OCEAN King India

Home of Quality Products



## Small Series DJC

### DJ-C Series Precise Electronic Balance

Model	Capacity Division	Pan Size (mm)	Display Method	Power
DJ102C	100g/0.01g	◆ 85	LCD Display	6*AAA 1.5V
DJ202C	200g/0.01g			
DJ302C	300g/0.01g			
DJ502C	500g/0.01g			
DJ602C	600g/0.01g			
DJ1002C	1000g/0.01g			

### Common Features:

- \* Large LCD display with auto backlight
- \* Auto-tare
- \* Selectable measure units:g.oz.kg.lb.tl.dwt
- \* Overload protection
- \* Scale dimension: 155\*130\*55 (mm)

#### Contact Address :

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

#### Sales and Support:

Tel : +91-161-5049059, 2222471  
Email: info@oceankingindia.com, amitsood24@gmail.com  
website: www.oceankingindia.com, www.goldtester.co