



# OCEAN King India

Home of Quality Products



**Desktop Kind 7200:** This high tech system is being designed and developed by our manufacturers in Superiority with global quality standards.

### **Technical Specification for Desktop Kind 7200 Gold Testing Machine**

- Analysis Elements: Au(Gold), Ag(Silver), Pt(Platinum), Pd(Palladium), Cu(Copper), Zn(Zinc), Ni(Nickel).
- Analysis Range: 0.2%.99%
- Test point dimension: 12mm
- Sample Case Size: 090 mm
- Spectrometer Size: L500 mm\*W450mm\*H350 mm
- Test Time: Confirm the elements in 13 seconds, operator could set the time for testing the elements indensity)
- Test Precision: +/-0.1% The Precision could keep 0.1 after several repeated test
- High Pressure Template 450Kv
- X-Ray Source: Mo Material X-ray light tube (Wind Cool, No radiation)
- Max. Power: 120W
- Weight: 35Kg
- Detector: Hermitic Proportional Counter with Micro processor
- Coated layer measurement: Coated layer thickness Range <30µm
- Eight fans of the circulatory system, Temperature control system can protect themselves.

#### **Contact Address :**

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

#### **Sales and Support:**

**Tel :** +91-161-5049059, 2222471

**Email:** info@oceankingindia.com, amitsood24@gmail.com  
website: www.oceankingindia.com, www.goldtester.co



**OCEAN King India**  
Home of Quality Products

**Software:**

- A Powerful, yet easy to use Software, with various features.
- Powerful templates created for easy reference and standards comparison
- Data transit between systems
- Data backup
- Hardware and Software lock is provided to prevent any kind of misuse by unauthorized person
- Operation and Administrator mode to limit and delimit access
- Easy to initialize by the user in one single step with Standard Silver Reference Material
- Operating system – Windows 2000, XP & Vista
- Spectrum Analysis with real time and pre-defined similar samples to analyze the differential results
- Reports can be generated on various parameters to suit the client
  
- Certificates of ornaments with Image can be generated

**Other Information:**

- **Pay Mode Terms:** L/C (Letter of Credit), TT Bank Transfer, Western union, Money gram
- **Port of Dispatch:** New Delhi, India
  
- **Packaging Details:** International Export Standard packaging.

**Contact Address :**

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

**Sales and Support:**

**Tel :** +91-161-5049059, 2222471

**Email:** info@oceankingindia.com, amitsood24@gmail.com  
website: www.oceankingindia.com, www.goldtester.co