



# OCEAN King India

Home of Quality Products



**Portable Kind 6308** is the New Portable model design for Customers who need Mobile Application. The outer cover has been changed to nowadays the popular anti-collision engineering plastics, Durable with Elegant appearance. This new model equipped High resolution digital camera. Portable model is strongly recommended to the Jewelry recycling, Gold purchase from Mining companies or traders and the field that need to carry the instrument to travel around. Test Gold Dust Ore Bars And Ingots

**Contact Address :**

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

**Sales and Support:**

**Tel :** +91-161-5049059, 2222471

**Email:** info@oceankingindia.com, amitsood24@gmail.com  
website: www.oceankingindia.com, www.goldtester.co



# OCEAN King India

Home of Quality Products

## Detailed Features:

Analysis Elements: Au(Gold), Ag (Silver), Pt (Platinum), Pd(Palladium), Cu(Copper), Zn(Zinc), Ni(Nickel)

Analysis Range:10% ~ 99.99%

Test point dimension:1 ~ 2mm

Sample Chamber Size: 0 ~ 90 mm

Spectrometer Size: L500 mm\*W400mm\*H200 mm

Test Time: Confirm the elements in 1 ~ 3 seconds, operator could set the time

Test Precision:+/-0.2% (The Precision could keep 0.2 % after several repeated test)

High Voltage:4 ~ 50Kv

X-Ray Source: Mo Material X-ray light tube (Wind Cool, No radiation)

Max. Power: 100W

Weight: 20Kg

Detector: Hermitic Proportional Counter with Micro processor

Coated layer measurement: Coated layer thickness Range <30µm

Work Temperature: 5-46°C

Software: XRF Version 3.0

Operation system: Windows2000/Me/XP

## Other Information:

- Pay Mode Terms: L/C (Letter of Credit),TT Bank Transfer, Western union, Money gram
- Port of Dispatch: New Delhi, India
- Packaging Details: International Export Standard packaging.

### Contact Address :

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

### Sales and Support:

Tel : +91-161-5049059, 2222471

Email: info@oceankingindia.com, amitsood24@gmail.com  
website: www.oceankingindia.com, www.goldtester.co