



**OCEAN** King India  
Home of Quality Products



**Ocean King Stand Alone 7800** is a Stand alone integrated XRF gold testing instrument with inbuilt Monitor shows you the result instantly.

## FEATURES

- Compact & Neat Design - Just One Power Cable & Its Ready to Work.
- Communication Back Port Like USB and RJ45 for. Stand Alone 7800 is ready to use with Wifi or Local Area Connection or VPN(Virtual Private Network)
- Beautiful appearance with Integrated PC & Monitor saves premium space in Jewellery Showrooms
- Easy to Use software means no confusion for users.
- Laser guided placement spot makes placing jewellery inside the chamber quick & easy.
- No more adjustments of stage or hovering in the test chamber to locate excitation

### Contact Address :

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

### Sales and Support:

**Tel :** +91-161-5049059, 2222471  
**Email:** info@oceankingindia.com, amitsood24@gmail.com  
website: www.oceankingindia.com, www.goldtester.co



# OCEAN King India

Home of Quality Products

- An important advantage of Bottom up X-Ray mounting.  
& Many more features, advantages & benefits apart from being Most Competitively Priced.

Specially Configured software & our technical knowhow helps the user get most accurate results. After sales service is the key, Ocean King has complete support to its customers and gives 100% backup across India, Europe, Canada, USA. Our engineers are Well trained & We stock good levels spare parts to support customers.

Professionally designed with Multi-node pre-start installations and the real- time monitoring system, the Stand Alone 7800 is an intelligent design with Modern Aesthetics and appeal. In order to ensure the user security, it's installed the automated anti-radiation device. Stand Alone 7800 gold testing instrument is ideal for showroom, manufacturing centers, wholesalers, Banks for gold loan, Export oriented jewellery manufacturing units, etc. to provide fast and accurate results.

#### **Detailed features:**

- Analysis Elements: Au, Ag, Pt, Pd, Rh, Ru, Cu, Zn, Ni, Cd, etc.
- Analysis Range: 1% 99.99%
- Test point dimension: 1 2mm
- Test Time: Confirm the elements in 30 to 200 second, operator could set the time for testing the elements intensity.)
- Coated layer measurement: Coated layer thickness Range <30µm
- Test Precision: +/-0.1% The Precision could keep 0.1 after several repeated test)
- High Voltage Generator: 450Kv
- X-Ray Source: Mo Material X-ray light tube (Wind Cool, No radiation)
- Detector: High Resolutions Hermetic Proportional Counter with Micro processor technology Auto thermal environ change cut off system, protects the system from Damage or Malfunctions, thus protecting the results derived on the system.
- Integrated PC with in-built monitor, designed compact with Small foot print.
- Spectrometer Size: L500 mm\*W500mm\*H400 mm
- Max. Power: 120W
- Weight: 50Kg
  
- Temp Control: Eight fans of the circulatory system, Temperature control system can protect themselves

#### **Contact Address :**

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

#### **Sales and Support:**

**Tel :** +91-161-5049059, 2222471

**Email:** info@oceankingindia.com, amitsood24@gmail.com

**website:** www.oceankingindia.com, www.goldtester.co



**OCEAN King India**  
Home of Quality Products

### **Software:**

#### **A Powerful, yet easy to use Software, with various features:**

- Powerful templates created for easy reference and standards comparison
  - Data transit between systems.
  - Integrated PC – saves space due to small foot print.
  - Hardware and Software lock is provided to prevent any kind of misuse by unauthorized person.
  - Operation and Administrator mode to limit and delimit access.
  - Easy to initialize by the user in one single step with Standard Silver Reference Material.
- 
- Operating system – Windows 2000, XP & Vista

### **Other Information**

- Minimum Order Quantity: 1 Set/Sets
  - Port of Dispatch: New Delhi, India
  - Delivery Time: 2 - 3 weeks
- 
- Packaging Details: Wooden crate packing for Safe & Secured delivery across India Standard Export Packing.

#### **Contact Address :**

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

#### **Sales and Support:**

**Tel :** +91-161-5049059, 2222471

**Email:** info@oceankingindia.com, amitsood24@gmail.com  
website: www.oceankingindia.com, www.goldtester.co