



OCEAN King India

Home of Quality Products



Stand Alone 9500

Dual display screen design, professional fashion appearance with multi-node pre-start installations and the real-time monitoring system; The cover locks and electronic locks are making day-to-day management more convenient;

EXF9500 have integrated industrial computers, hence, no need to connect with an external computer; Its also installed the automated anti-radiation device in order to ensure the user security. Electric devices for opening cover makes the test easier. Ultra-high resolving power can easily settle any complexity requirements for non-destructive analysis the in precious metals industry. Among the superlative spectrometers in precious metals industry. EXF9500 is the best competitive one.

Common Features:

- * **Analysis Elements:** Au (Gold), Ag (Silver), Pt (Platinum), Pd (Palladium), Rh (Rhodium), Ru (Ruthenium), Cu (Copper), Zn (Zinc), Ni (Nickel), Cd Cd (Cadmium), In (Indium), Os (Osmium), Ir (Iridium)
- * **Analysis Range:** 0.01%99.99%
- * **Test point dimension:** 12mm
- * **Sample Chamber Size:** 090 mm
- * **Spectrometer Size:** L500 mm*W500mm*H450 mm
- * **Test Time:** Confirm the elements in 13 seconds, operator could set the time for testing the elements in density.
- * **Test Precision:** +/-0.01% (The Precision could keep 0.01 after several repeated test)

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471
Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co



OCEAN King India

Home of Quality Products

- * **High Pressure Template:** 450Kv
- * **X-Ray Source:** Mo Material X-ray light tube(Wind Cool, No radiation)
- * **Max. Power:** 144W
- * **Weight:** 60Kg
- * **Detector:** Si-Pin (Cooling: Semiconductor)
- * Eight fans of the circulatory system, Temperature control system can protect themselves

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471

Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co