




OCEAN King India

Home of Quality Products



Specially Configured software & our technical knowhow helps the user get most accurate results. After sales service is the key, OceanKing has complete support to its customers and gives 100% backup across India. Our engineers are well trained & We stock good levels of spare parts to support our customers.

Ocean King Stand Alone 8000S is excellent in design & delivers powerful performance. Its an excellent product offering highest Price to Performance Ratio in the industry. The New Design is based on Aero Space Design concept and also its sleek shape makes it a compelling instrument for Every Jewellery Showroom Or Banks. In order to ensure the user safety, it's installed the automated anti-radiation device. Stand Alone 8000S gold testing instrument is ideal for showroom, manufacturing centers, wholesalers, Banks for gold loan, Export oriented jewellery manufacturing units, etc. to provide fast and accurate results.

Detailed Features:

- Analysis Elements: Au(Gold), Ag(Silver), Pt(Platinum), Pd(Paladium), Rh(Rohdium), Ru(Ruthenium), Cu(Copper), Zn(Zinc), Ni(Nickel), Cd(Cadmium), etc.

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471
Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co



OCEAN King India

Home of Quality Products

- Analysis Range: 0.1% 99.99%
- Test point dimension: 1.2mm
- Test Time: Confirm the elements in 30~200 second, operator could set the time for testing the elements intensity.
- Coated layer measurement: Coated layer thickness Range <math><30\mu\text{m}</math>
- Test Precision: $\pm 0.1\%$ The Precision could keep 0.1 after several repeated test)
- X-Ray Source: Mo Material X-ray light tube (Wind Cool, No radiation)
- Detector: High Resolutions Hermetic Proportional Counter with Micro processor technology Auto thermal environ change cut off system, protects the system from Damage or Malfunctions, thus protecting the results derived on the system.
- Integrated PC with in-built monitor, designed compact with Small foot print.

- Temp Control: Eight fans of the circulatory system, Temperature control system can protect themselves

Software:

A Powerful, yet easy to use Software, with various features:

- Powerful templates created for easy reference and standards comparison
- Data transit between systems.
- Integrated PC – saves space due to small foot print.
- Hardware and Software lock is provided to prevent any kind of misuse by unauthorized person.
- Operation and Administrator mode to limit and delimit access.
- Easy to initialize by the user in one single step with Standard Silver Reference Material.
- Operating system – Windows 2000, XP & Vista

- Standard Holders for Round Samples Like Ring and Bangles

Other Information

- Pay Mode Terms: T/T (Bank Transfer)
- Minimum Order Quantity: 1 Set/Sets
- Port of Dispatch: New Delhi India

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471

Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co



OCEAN King India

Home of Quality Products

- Delivery Time: 2 Weeks
- Packaging Details: Wooden crate packing for Safe & Secured delivery across India and Standard Export Packing

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471

Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co