



OCEAN King India

Home of Quality Products



Ocean King Gold Testing Machine Desktop Kind 9500 is based on Si-Pin Diode Detection technology. Ocean King Desktop Kind 9500 is capable of detecting all metals. Due to its high resolution detection capability it detects any metal in lowest PPM concentration. Thus low content elements are easily analyzed and detected for its presence in the Alloy. This is essential function for Assay Labs & Testing Centers which provide certificates for tested sample. It's also installed with automated anti-radiation device in order to ensure the user safety.

Touch Screen LCD control panel is fitted in the instrument to provide realtime system status to the user, this usually helps the user to control the essential instrument parameters such as Temperature and power supply. With user friendly and comprehensive analysis software design. The user is at a great convenience in deriving accurate results in minimum time and least efforts. Result display in both percentage and karat. Automatically repeat testing multiple times and get average results. Multiple test report format to meet different customer's printing requirement. Desktop Kind 9500 is best for Hallmarking centers, Assaying Centers, Testing Centres, etc. Where additional elements required to be tested, where certification of sample is the task. Desktop Kind 9500 is the answer.

**Analysis elements: AU(Gold) AG(Silver) PT(Platinum) PD(Palladium)
RH(Rohdium) RU(Ruthenium) CU(Copper) ZN (Zinc) NI(Nickel) CD(Cadmium) in
OS IR(Iradium) W (from Na to U, around 79 elements).**

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471
Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co



OCEAN King India

Home of Quality Products

- **Analysis range:** 0. 01% 99. 99%
- **Test time:** Typical 100 seconds (operator could set the time according to situation)
- **Test point dimension:** 1~2mm
- **Test precision:** +/-0. 01%-0. 1%
- **Instrument current:** 120w
- **High voltage template:** 50kv
- **X-ray source:** mo material x-ray light tube (wind cool, no radiation)
- **Detector:** si-pin (cooling: semiconductor)
- **Coated layer measurement:** coated layer thickness range <30µm
- **Computer :** Externally connected Both PC / Laptop.USB Port Communication Port
- **Operation system:** windows 7

- **Net weight:** 30kg

Applications:

- Assayers
- Hallmark labs
- Testing centers

- Refiners & high precise analysis applications.

Other Information

- **Pay Mode Terms:** Bank Transfer (TT) Other
- **Minimum Order Quantity:** 1 Set/Sets
- **Port of Dispatch:** New Delhi , India

- **Delivery Time:** 2 Weeks

Packaging Details: Standard Wooden Crate Packing.

Contact Address :

Ocean King International
6213/B, St. No.-4, Hargobind Nagar,
Ludhiana, India 141008.

Sales and Support:

Tel : +91-161-5049059, 2222471
Email: info@oceankingindia.com, amitsood24@gmail.com
website: www.oceankingindia.com, www.goldtester.co