



# OCEAN King India

Home of Quality Products



**Ocean King Stand Alone 9600U** Integrated PC based with Multitouch Touch Screen Interface. SI Pin Diode Silicon Pin Detector It's a stand alone integrated XRF gold testing instrument with inbuilt Muti Touch Screen Laptop shows you the result instantly. With smaller foot print and due to integrated System, the unit saves space time.

## FEATURES

- Compact Design without Messy Cables
- Beautiful appearance with Integrated PC & Monitor saves space
- Easy to Use GUI software means no confusion for users.XRF 6.0
- Laser guided placement spot makes placing jewellery inside the chamber quick & easy.
- No more adjustments of stage or hovering in the test chamber to locate excitation - An important advantage of Bottom up X-Ray mounting.

& Many more features, advantages & benefits apart from being Most Competitively Priced Gold Karat Analyzer In India.

### Contact Address :

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

### Sales and Support:

Tel : +91-161-5049059, 2222471  
Email: info@oceankingindia.com, amitsood24@gmail.com  
website: www.oceankingindia.com, www.goldtester.co



# OCEAN King India

Home of Quality Products

Specially Configured software & our technical knowhow helps the user get most accurate results. After sales service is the key, Ocean King has complete support to its customers and gives 100% backup across India. Our engineers are Well trained & We stock good levels of spare parts to support our customers.

Ocean King Stand Alone 9600U is excellent in design & delivers powerful performance. Its an excellent product offering highest Price to Performance Ratio in the industry. The New Design is based on Aero Space Design concept and also its sleek shape makes it a compelling instrument for

Every Jewellery Showroom Or Banks. In order to ensure the user safety, it's installed the automated anti-radiation device. Stand Alone 9600U gold testing instrument is ideal for showroom, manufacturing centers, wholesalers, Banks for gold loan, Export oriented jewellery manufacturing units, etc. to provide fast and accurate results.

**Analysis elements: AU(Gold) AG(Silver) PT(Platinum) PD(Palladium) RH(Rohdium RU(Ruthenium CU(Copper) ZN (Zinc) NI(Nickel) CD(Cadmium) in OS IR(Iradium) W (from Na to U, around 79 elements).**

- **Analysis range:** 0. 01% 99. 99%
- **Test time:** Typical 100 seconds (operator could set the time according to situation)
- **Test point dimension:** 1~2mm
- **Test precision:** +/-0. 01%-0. 1%
- **Instrument current:** 120w
- **High voltage template:** 50kv
- **X-ray source:** mo material x-ray light tube (wind cool, no radiation)
- **Detector:** si-pin (cooling: semiconductor)
- **Coated layer measurement:** coated layer thickness range <30µm
- **Computer :** Internal Touch Screen Laptop. Inbuild Thermal Printer
- **Operation system:** windows 7
  
- **Net weight:** 30kg

**Contact Address :**

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

**Sales and Support:**

**Tel :** +91-161-5049059, 2222471

**Email:** info@oceankingindia.com, amitsood24@gmail.com  
website: www.oceankingindia.com, www.goldtester.co



**OCEAN King India**  
Home of Quality Products

**Applications:**

- Assayers
- Hallmark labs
- Testing centers
  
- Refiners & high precise analysis applications.

**Other Information**

- **Pay Mode Terms:** TT Wire Transfe, Other
- **Minimum Order Quantity:** 1 Set/Sets
- **Port of Dispatch:** New Delhi, India
- **Delivery Time:** 2 Weeks
  
- **Packaging Details:** Standard Wooden Crate Packing

**Contact Address :**

Ocean King International  
6213/B, St. No.-4, Hargobind Nagar,  
Ludhiana, India 141008.

**Sales and Support:**

**Tel :** +91-161-5049059, 2222471

**Email:** info@oceankingindia.com, amitsood24@gmail.com  
website: www.oceankingindia.com, www.goldtester.co